

Instrumentation class, SMS 598, Fall 2012
Emmanuel Boss

Group I, attenuation lab:

In this lab you will build an attenuation meter as well as use a commercial one. You will use it to measure the attenuation of food coloring (which is the same as its absorption as it is a dissolved substance that negligibly scatters).

Follow the instructions in the associated worksheet to construct a transmissometer with which you will determine the specific attenuation of a drop of food coloring. In parallel you will be using a commercial beam attenuation meter (WETLabs C-star). If time permits run filtered and total seawater sample and compare with ac-S value for the same wavelength.

Appendix: The Beer-Lambert-Bouguer law:

The law is an empirical law relating the attenuation of light to the properties of the material the light is traveling through. It states that when we shine a collimated light beam through a liquid and measure the light getting to a receiver at distance r , the intensity of the light getting to the receiver ($I(r)$) is a function of the light leaving the source (I_0), the concentration of the absorbing (or attenuating, for material that scatter as well as absorb) substance, $[c]$ the mass (or mol) specific absorption (or attenuation) by the substance, c^* , and the path-length, r :

$$I(r)=I_0*\exp(-c^*[c]r). \quad (1)$$

Beer's Law often refers to absorption by solutes while Lambert's (and Bouguer) law refers to attenuation by suspension of particles in a liquid (e.g. Shifrin, 1988, <http://scienceworld.wolfram.com/physics/LambertsLaw.html>, and <http://scienceworld.wolfram.com/physics/BeersLaw.html>).

The nondimensional quantity: $2.303 c^*[c]r$ is often referred to as the 'Absorbance' when one deals with purely absorbing substances or 'Attenuance' when one deals with attenuating (absorbing+scattering) material. $2.303=\log_e(10)$.

The law fails at high concentration when either:

1. Error due to detector sensitivity become significant (absorbing materials).
2. Multiple scattering causes deviations from Equ. 1 (scattering materials).

History (From <http://www.worldhistory.com/wiki/B/Beer-Lambert-law.htm>):

Beer's law was independently discovered (in various forms) by [Pierre Bouguer](#) (Mathematician and hydrographer, France) in [1729](#), [Johann Heinrich Lambert](#) (Physicist and astronomer, Germany, other contribution to optics abound, in particular the law about reflection from surfaces) in [1760](#) and [August Beer](#) (Mathematician, chemist and physicist, Germany) in [1852](#).

While this law was discovered empirically it turns out it can be derived from theory. The statistics of interaction of light with randomly distributed molecules/particles give rise to an exponential attenuation within a substance (Kostinski, 2001).

Kostinski, A. B., 2001. On the extinction of radiation by a homogeneous but spatially correlated random medium. JOSA A, Vol. 18, Issue 8, pp. 1929-1933.

Shifrin, K. 1988. Physical Optics of Ocean Water. AIP. 285pp.