

Recap laboratory exercise, July 13, 2015

**AC-meters: measure Milli-Q, filtered DRE and total DRE**

For absorption tube, plot:

raw data

Milli-Q corrected dissolved and total absorption

T corrected dissolved and total absorption AND dissolved and total attenuation

S corrected dissolved and total absorption AND dissolved and total attenuation

Total absorption minus dissolved absorption (both TS corrected).

Spectra for each scatter-corrected particulate absorption (different versions, one graph per version).

**Spectrophotometer:**

Wavelengths: 350 – 800 nm, slow scan speed, zero against air.

Pad samples:

Filter 100 mL seawater for blank – do one filter

Filter 400 mL for sample – do replicates

Seawater blank:

Plot raw blank pad spectrum

Plot blank pad spectrum with null value in IR subtracted

DRE particulate spectrum:

Scan pads with DRE water

Plot raw a\_part pad spectrum with null value in IR subtracted

Plot a\_part scan with blank spectrum subtracted

Methanol extract:

Pour hot methanol on blank and one DRE filter

Allow to sit for 10 min; add more hot methanol if necessary

Scan blank and DRE pad, and plot each step as above.

Derive a-phyt as the difference between fully corrected spectra for pads before and after methanol treatment.

Compare a\_part( $\lambda$ ) between ac-meter and spectrophotometer.