

SOCCOM P18 HPLC

Nils Haëntjens¹, Emmanuel Boss¹, and Lynne Talley²
¹University of Maine, ²Scripps Institution of Oceanography
October 18, 2017

Sample collection

Near-surface samples from SOCCOM CTD stations were taken for HPLC analysis. 1-2 L of sample was filtered in the dark through glass fiber filter (GF/F) having a diameter of 25 mm. Filters were immediately stored in aluminium foil packages in a dewar of liquid nitrogen (-80 °C). Samples were shipped and analyzed at NASA GSFC.

More information on the cruise are available at:

<https://socom.princeton.edu/content/shipboard-data-reports>

Analysis method

The HPLC analysis method can be cited as Van Heukelem and Thomas (2001), further described in Hooker et al. (2005). For a more detailed description, please see below; contact Crystal Thomas (crystal.s.thomas@nasa.gov) for a tailored description.

The HPLC used for pigment analysis is an Agilent RR1200 with a programmable autoinjector (900 ul syringe head), refrigerated autosampler compartment, thermostatted column compartment, quaternary pump with in-line vacuum degasser, and photo-diode array detector with deuterium and tungsten lamps. The HPLC is controlled by Agilent Chemstation software.

The 4.6 x 150 mm HPLC Eclipse XDB column (Agilent Technologies, Palo Alto, CA) is filled with a C8 stationary phase (3.5 um stationary phase); the mobile phase consists of a linear gradient from 5-95% solvent B over 27 minutes, for which solvent A is 70 parts methanol, 30 parts 28 mM tetrabutylammonium acetate (pH 6.5) and solvent B is methanol. The column temperature is 60 °C and the photo diode array detector is set to plot chromatograms at 450, 665, and 222 nm to acquire visible absorbance spectra between 350 and 750 nm.

Vitamin E acetate is used as the internal standard (ISTD) for determining extraction volumes. Its absorbance is monitored at 222 nm; it has negligible absorbance at 450 nm and none at 665 nm. Therefore, it does not interfere at wavelengths used to quantify pigments and can be used in very high concentrations with S:N ratios much higher than are possible with pigments. The high signal:noise ratio contributes to excellent analysis precision, for which injection repeatability averages 0.6%. It is stable under conditions of extraction and analysis.

Calibration is performed with individual pigment standards, whose concentrations have been determined spectrophotometrically using absorption coefficients in common with those used by most other laboratories (Hooker et al. 2005) and the commercial vendor, DHI Water and Environment (Horsholm, Denmark). Standards are either purchased from DHI (in solution with concentrations provided) or purchased in solid form and suspended in solvent at GSFC. Thirty-six peaks are individually quantified by HPLC, from which 26 pigments are reported (some pigments contain individual components that are summed and reported as one pigment).

Abbreviations

Primary Pigments

| | |
|----------------|----------------------------|
| Allo | alloxanthin |
| alpha-beta-Car | carotenes |
| But-fuco | 19'-butanoyloxyfucoxanthin |
| Diadino | diadinoxanthin |
| Diato | diatoxanthin |
| Fuco | fucoxanthin |
| Hex-fuco | 19'-hexanoyloxyfucoxanthin |
| Perid | Peridinin |
| Tot_Ch_l_a | total chlorophyll a |
| Tot_Ch_l_b | total chlorophyll b |
| Tot_Ch_l_c | total chlorophyll c |
| Zea | Zeaxanthin |

Secondary Pigments

| | |
|-----------|---|
| Chl_c3 | Chlorophyll c3 |
| Chlide_a | chlorophyllide a |
| DV_Ch_l_a | divinyl chlorophyll a |
| DV_Ch_l_b | divinyl chlorophyll b |
| MV_Ch_l_a | monovinyl chlorophyll a |
| MV_Ch_l_b | monovinyl chlorophyll b |
| | Chlorophyll c2 + chlorophyll c1 + MGDVP |
| | Mg-2,4-divnyl pheoporphyrin a5 monomethyl ester |

Tertiary Pigments

| | |
|----------|----------------------|
| Lut | Lutein |
| Neo | Neoxanthin |
| Phide_a | total pheophorbide a |
| Phytin_a | total pheophytin a |
| Pras | Prasincoxanthin |
| Viola | Violaxanthin |

Ancillary Pigment

| | |
|------|---------------------|
| Gyro | Gyroxanthin diester |
|------|---------------------|

Other abbreviations

| | | |
|----------|---|---|
| DP | total diagnostic pigments | PSC + allo + zea + Tot_Ch_l_b |
| PPC | photoprotective carotenoids | allo + diadino + diato + zea + alpha-beta-car |
| PPC_TCar | ratio of photoprotective carotenoids to total carotenoids | [PPC]/[Tcar] |
| PPC_TPg | ratio of photoprotective carotenoids to total pigments | [PPC]/[Tpg] |
| PSC | photosynthetic carotenoids | but-fuco + fuco + hex-fuco + perid |
| PSC_TCar | ratio of photosynthetic carotenoids to total carotenoids | [PSC]/[TCar] |
| PSP | photosynthetic pigments | PSC + TChl |
| PSP_TPg | ratio of photosynthetic pigments to total pigments | [PSP]/[TPg] |

| | | |
|------------|--|--------------------------------------|
| TAcc | total accessory pigments | PPC + PSC + Tot_Ch_l_b + Tot_Ch_l_c |
| TAcc_TChla | ratio of total accessory pigments to total chlorophyll a | [Tacc]/[Tchla] |
| TCar | total carotenoids | PPC + PSC |
| TChl | total chlorophylls | Tot_Ch_l_a + Tot_Ch_l_b + Tot_Ch_l_c |
| TChl_TCar | ratio of total chlorophyll to total carotenoids | [TChl]/[TCaro] |
| TChla_Tpg | ratio of total chlorophyll a to total pigments | [TChla]/[TPg] |
| TPg | total pigments | TAcc + Tot_Ch_l_a |