SBE 50

Digital Oceanographic Pressure Sensor



Note: NEW ADDRESS as of January 18, 2010

CE

<u>User's Manual</u>

Sea-Bird Electronics, Inc. **13431 NE 20th Street Bellevue, Washington 98005 USA** Telephone: 425/643-9866 Fax: 425/643-9954 E-mail: seabird@seabird.com Website: www.seabird.com

Manual Version #009, 01/29/10 Firmware Version 1.0d and later

Limited Liability Statement

Extreme care should be exercised when using or servicing this equipment. It should be used or serviced only by personnel with knowledge of and training in the use and maintenance of oceanographic electronic equipment.

SEA-BIRD ELECTRONICS, INC. disclaims all product liability risks arising from the use or servicing of this system. SEA-BIRD ELECTRONICS, INC. has no way of controlling the use of this equipment or of choosing the personnel to operate it, and therefore cannot take steps to comply with laws pertaining to product liability, including laws which impose a duty to warn the user of any dangers involved in operating this equipment. Therefore, acceptance of this system by the customer shall be conclusively deemed to include a covenant by the customer to defend, indemnify, and hold SEA-BIRD ELECTRONICS, INC. harmless from all product liability claims arising from the use or servicing of this system.

Table of Contents

Section 1: Introduction	4
About this Manual	4
Quick Start	
Unpacking SBE 50	5
Section 2: Description of SBE 50	6
System Description	6
Specifications	
Dimensions and End Cap Connector	9
Section 3: Preparing SBE 50 for Deployment	10
Software Installation	
Power and Communications Test	
Test Setup	
Test	11
Section 4: Deploying and Operating SBE 50	14
Sampling Modes	
Polled Sampling	
Autonomous Sampling	15
Baud Rate, Cable Length, Power, and Data Transmission Rate	
Baud Rate, Cable Length, and Data Transmission Rate	
Power and Cable Length	
Optional Voltage Output Interface Command Descriptions	
Data Output Formats	
Deployment	
Recovery	
Section 5: Routine Maintenance and Calibration	26
Corrosion Precautions	
Connector Mating and Maintenance	
Pressure Sensor Maintenance	
Sensor Calibration	
Section 6: Troubleshooting	29
Problem 1: Unable to Communicate with SBE 50	
Problem 2: Unreasonable Data	
Glossary	
Appendix I: Functional Description	
Sensor	
Sensor Interface	
Appendix II: Electronics Disassembly/Reassembly	
Appendix III: Command Summary	
Appendix IV: Replacement Parts	
Index	

Section 1: Introduction

This section includes a Quick Start procedure, and photos of a standard SBE 50 shipment.

About this Manual

This manual is to be used with the SBE 50 Digital Oceanographic Pressure Sensor. It is organized to guide the user from installation through operation and data collection. We've included detailed specifications, command descriptions, maintenance and calibration information, and helpful notes throughout the manual.

Sea-Bird welcomes suggestions for new features and enhancements of our products and/or documentation. Please contact us with any comments or suggestions (seabird@seabird.com or 425-643-9866). Our business hours are Monday through Friday, 0800 to 1700 Pacific Standard Time (1600 to 0100 Universal Time) in winter and 0800 to 1700 Pacific Daylight Time (1500 to 0000 Universal Time) the rest of the year.

Quick Start

Follow these steps to get a Quick Start using the SBE 50. The manual provides step-by-step details for performing each task:

- 1. Perform pre-check procedures (see *Section 3: Preparing SBE 50 for Deployment*):
 - A. Test power and communications.
- 2. Deploy the SBE 50 (see Section 4: Deploying and Operating SBE 50):
 - A. Establish setup and operating parameters.
 - B. Check status (DS) and calibration coefficients (DCal) to verify setup.
 - C. Install I/O cable connector and locking sleeve.
 - D. Deploy SBE 50, using customer-supplied hardware.

Unpacking SBE 50

Shown below is a typical SBE 50 shipment.





I/O Cable (for standard 4-pin connector)



SBE 50 User Manual



Software, and Electronic Copies of Software Manuals and User Manual

Section 2: Description of SBE 50

This section describes the functions and features of the SBE 50, including specifications and dimensions.

System Description

The SBE 50 is a high-accuracy, high-resolution, high sample rate (16 Hz sampling) pressure monitor, designed to be incorporated into oceanographic sensing systems at depths up to 7000 meters (23,000 feet).

The SBE 50's pressure sensor is offered in eight full-scale ranges from 20 to 7000 decibars. Developed by Druck, Inc., the pressure sensor has a superior new design that is entirely different from conventional *silicon* types in which the deflection of a metallic diaphragm is detected by epoxy-bonded silicon strain gauges. The Druck sensor employs a micro-machined *silicon diaphragm* into which the strain elements are implanted using semiconductor fabrication techniques. Unlike metal diaphragms, silicon's crystal structure is perfectly elastic, so the sensor is essentially free of pressure hysteresis. Compensation of the temperature influence on pressure offset and scale is performed by the SBE 50's CPU.

The SBE 50 is an easy-to-use, light, compact instrument ruggedly made of titanium and other low-maintenance (plastic) materials; it is well suited to even the smallest vehicle. Commands sent to the SBE 50 provide status display, data acquisition setup, and diagnostic tests. Its data is telemetered by its internal, 3-wire, RS-232C interface. The SBE 50 has no memory or internal batteries.

The SBE 50 samples in two modes:

- Autonomous sampling –The SBE 50 runs continuously and samples at 16 Hz (16 scans per second). It can be set to average up to 255 samples, and transmits the averaged data. The SBE 50 can be programmed to begin autonomous sampling when power is applied or on command.
- **Polled sampling** On command, the SBE 50 takes one sample and transmits the data.


The SBE 50 outputs data in one of the following formats:

- Raw pressure and pressure temperature values.
- Data in ASCII engineering units (pressure in decibars or psia, or depth in meters or feet) pressure is calculated using calibration coefficients stored in EEPROM; depth is calculated from pressure, based on user-selection of fresh or salt water and user-input latitude (for salt water only).
- Pressure and scan number in hexadecimal this format is used when the SBE 50 interfaces with the Auto Fire Module (AFM) for controlling an SBE 32 Carousel Water Sampler, or when the SBE 50 interfaces with an SBE 55 ECO Water Sampler.

Options include:

- Voltage Output interface, which provides 0-5 volt output proportional to the pressure or depth, *in addition to* the RS-232C output.
- An MCBH Micro (wet-pluggable) connector in place of the standard glass-reinforced epoxy connector.

The SBE 50 can be used with the following Sea-Bird equipment:




SBE 32 Carousel Water Sampler and Auto Fire Module (AFM) -The SBE 32 is a 12-, 24, or 36-bottle water sampler, available with bottle sizes up to 30 liters. The AFM, mounted in a bottle position on the SBE 32, allows the SBE 32 to operate autonomously on non-conducting cables. The AFM supplies the power, logic, and control commands to operate the SBE 32, and supplies power for the SBE 50. The AFM monitors the pressure data transmitted by the SBE 50 in real-time, closing bottles at predefined pressures (depths) or whenever the system is stationary for a specified period of time. Bottle number, firing confirmation, and five scans of SBE 50 data are recorded in the AFM memory for each bottle fired. See the AFM manual for details.

SBE 55 ECO Water Sampler – The ECO is a 3- or 6-bottle water sampler with 4-liter bottles. The ECO has integrated electronics that function similar to the AFM. The ECO supplies power for the SBE 50 and monitors the pressure data transmitted by the SBE 50 in real-time, closing bottles at predefined pressures (depths) or whenever the system is stationary for a specified period of time. Bottle number, firing confirmation, and five scans of SBE 50 data are recorded in the ECO memory for each bottle fired. See the ECO manual for details.

Note:

Seasave (our real-time data acquisition software) does not support the SBE 50 used with the SBE 36/PDIM. Real-time data can be viewed in SEATERM (terminal program). Surface PAR data (if applicable) is output as raw data. SBE 36 CTD Deck Unit and Power Data Interface Module (PDIM) – These items provide power and real-time data handling capability over single-conductor sea cables. The PDIM is a small pressure housing that is mounted with the SBE 50. It provides +15 VDC power to the SBE 50 and interfaces two-way RS-232 communications from the SBE 50 to the telemetry used on the sea cable. See the SBE 36/PDIM manual for operating details.



Notes:

- Help files provide detailed information on the software.
- See the AFM manual for details on the use of SeatermAF with an SBE 50 and AFM.
- See the SBE 55 ECO manual for details on the use of SeatermAF with an SBE 50 and ECO Water Sampler.
- Sea-Bird supplies the current version of our software when you purchase an instrument. As software revisions occur, we post the revised software on our FTP site. See our website (www.seabird.com) for the latest software version number, a description of the software changes, and instructions for downloading the software from the FTP site.

The SBE 50 is supplied with a powerful 2000/XP software package, SEASOFT V2, which includes:

- **SEATERM** terminal program for communication and data acquisition. SEATERM can send commands to the SBE 50 to provide status display, data acquisition setup, data display and capture, and diagnostic tests.
- SeatermAF terminal program for communication with the SBE 50 when used with an SBE 32 Carousel Water Sampler and AFM, or with an SBE 55 ECO Water Sampler.

Note that SEATERM and SeatermAF do not process the data.

Specifications


Note:

Pressure ranges are expressed in meters of deployment depth capability.

Measurement	0 to full scale range:		
Range	20 / 100 / 350 / 600 / 1000 / 2000/ 3500 / 7000 meters		
Initial Accuracy	0.1% of full scale range		
Typical Stability (per year)	0.05% of full scale range		
Resolution	0.002% of full scale range (D/A resolution for optional 0-5 volt output: 16 bit)		
Sensor Calibration	Ambient pressure to full scale range in 5 steps		
External Input Power	8-30 VDC		
_	Input (Volts) Current (milliamps)		
Power	8 51		
Requirements	10 43		
(with voltage	15 29		
output option	20 22		
installed)	25 19 30 21		
Housing	Titanium pressure case rated at 7000 meters		
	(23,000 feet)		
Weight	In air: 0.7 kg (1.5 lbs)		
theight	In water: 0.4 kg (0.9 lbs)		

Dimensions and End Cap Connector

Dimensions in millimeters (inches)



6 8-30 VDC power in

Section 3: Preparing SBE 50 for Deployment

This section describes installing software and testing power and communications.

Software Installation

Notes:

- Help files provide detailed information on the software.
- It is possible to use the SBE 50 without the SEATERM terminal program by sending direct commands from a dumb terminal or terminal emulator, such as Windows HyperTerminal.
- Sea-Bird supplies the current version of our software when you purchase an instrument. As software revisions occur, we post the revised software on our FTP site. See our website (www.seabird.com) for the latest software version number, a description of the software changes, and instructions for downloading the software from the FTP site.


Sea-Bird recommends the following minimum system requirements for installing the software: Windows 2000 or later, 500 MHz processor, 256 MB RAM, and 90 MB free disk space for installation. Although SEASOFT V2 was designed to work with a PC running Win 2000/XP; extensive testing has not shown any compatibility problems when using the software with a PC running Windows Vista.

If not already installed, install SEATERM and other Sea-Bird software programs on your computer using the supplied software CD:

- 1. Insert the CD in your CD drive.
- Install software: Double click on SeasoftV2_date.exe (date is the date that version of the software was created). Follow the dialog box directions to install the software. The installation program allows you to install the desired components. Install all the components, or just install SEATERM.

The default location for the software is c:\Program Files\Sea-Bird. Within that folder is a sub-directory for each program.

Power and Communications Test





The power and communications test will verify that the system works, prior to deployment.

Test Setup

- 1. By hand, unscrew the locking sleeve from the SBE 50's bulkhead connector. If you must use a wrench or pliers, be careful not to loosen the bulkhead connector instead of the locking sleeve.
- 2. Remove the dummy plug from the SBE 50's I/O bulkhead connector by pulling the plug firmly away from the connector.
- Standard Connector Install the I/O cable connector, aligning the raised bump on the side of the connector with the large pin (pin 1 common) on the SBE 50. OR
 MCBL Connector Install the I/O cable connector aligning the give

MCBH Connector – Install the I/O cable connector, aligning the pins.

- 4. Connect the I/O cable connector to your computer's serial port.
- 5. Connect the I/O cable connector's red (+) and black (-) wires to a power supply (8-30 VDC).


Test

1. Double click on SeaTerm.exe. If this is the first time the program is used, the setup dialog box may appear:



Select the instrument type (*SBE 50*) and the computer COM port for communication with the SBE 50. Click OK.

2. The main screen looks like this:



Note:

Note:

See SEATERM's help files for

detailed information.

There is at least one way, and as many as three ways, to enter a command:

- Manually type a command in Command/Data Echo Area
- Use a menu to automatically generate a command
- Use a Toolbar button to automatically generate a command

Note:

Once the system is configured and connected (Steps 3 through 5 below), to update the Status bar:

- on the Toolbar, click Status; or
- from the Utilities menu, select

Instrument Status. SEATERM sends the status command, which displays in the Command/Data Echo Area, and updates the Status bar.


- Menus Contains tasks and frequently executed instrument commands.
- Toolbar Contains buttons for frequently executed tasks and instrument commands. All tasks and commands accessed through the Toolbar are also available in the Menus. To display or hide the Toolbar, select View Toolbar in the View menu. Grayed out Toolbar buttons are not applicable.
- Command/Data Echo Area Echoes a command executed using a Menu or Toolbar button, as well as the instrument's response. Additionally, a command can be manually typed in this area, from the available commands for the instrument. Note that the instrument must be *connected* for it to respond to a command (click Connect on the Toolbar to connect the instrument).
- Status bar Provides status information. To display or hide the Status bar, select View Status bar in the View menu.

Following are the Toolbar buttons applicable to the SBE 50:

Toolbar Button	Description	Equivalent Command*
Connect	Re-establish communications with SBE 50.	(press Enter
Connect	Computer responds with S> prompt.	key)
Status	Display instrument setup.	DS
Coefficients	Display calibration coefficients.	DCal
Capture	Capture instrument responses on screen to file. File has .cap extension. As SBE 50 has no internal memory, you must capture before sampling begins to save data for future review and processing. Press Capture again to turn off capture. Capture status displays in Status bar.	
Disconnect	Free computer COM port used to communicate with SBE 50. COM port can then be used by another program. Note that SBE 50 must be connected to COM port for data to be obtained, since it has no internal memory.	_

*See Command Descriptions in Section 4: Deploying and Operating SBE 50.

3. In the Configure menu, select SBE 50. The dialog box looks like this:



Make the selections in the Configuration Options dialog box:

- COMM Port: COM 1 through COM 10, as applicable
- Baud Rate: 9600 (documented on Configuration Sheet) •
- Data Bits: 8

Note:

the Configuration Options settings to the SeaTerm.ini file in your Windows

directory. SeaTerm.ini contains the

SEATERM and select the desired

instrument (SBE 37, 50, etc.) in the

Configure menu, the Configuration Options dialog box shows the last saved settings for that instrument.

last saved settings for each

instrument. When you open

- Parity: None •
- Mode: RS-232 (Full Duplex)
- Click OK to save the settings.

- 4. In the Communications menu, select *Options / Cycle baud when connecting*.
- Click Connect on the Toolbar. SEATERM tries to connect to the SBE 50 at the baud set in Step 3. If it cannot, it cycles through all other possible baud rates to try to connect. When it connects, the display looks like this: S>

This shows that correct communications between the computer and the SBE 50 has been established.

If the system does not respond with the S> prompt:

- Click Connect again or press the Enter key twice.
- Verify the correct instrument was selected in the Configure menu and the settings were entered correctly in the Configuration Options dialog box. Note that the baud rate is documented on the Configuration Sheet.
- Check cabling between the computer and SBE 50.
- 6. Display SBE 50 status information by clicking Status on the Toolbar. The display looks like this:

```
SBE50 V 1.0b SERIAL NO. 0011
Number of scans to average = 1
Start sampling on power up = no
Output format = depth, salt, meters, 2 decimals
Voltage output enabled
```

7. Command the SBE 50 to take a sample by typing **TS** and pressing the Enter key. The display looks like this (with output formatted for depth for salt water in meters, 2 digits after the decimal point):

0.06

where 0.06 = depth in meters This number should be reasonable; i.e., corresponding to approximately atmospheric pressure (0 decibars, 14.7 psia, or 0 depth).

8. Send commands to set up the instrument as desired, after reviewing the information in *Section 4: Deploying and Operating SBE 50* on sampling modes and commands.

The SBE 50 is ready for deployment.

Section 4: Deploying and Operating SBE 50

Note:

See SEATERM's help files for detailed information.

This section provides instructions for deploying the SBE 50. It includes:

- system operation
- example sets of operation commands
- cable length limitations
- setup of optional voltage output interface
- detailed command descriptions
- data output formats
- deployment
- recovery

Sampling Modes

The SBE 50 has two basic sampling modes for obtaining data: Polled Sampling and Autonomous Sampling.

Shown below are descriptions and examples of the two basic sampling modes. Note that the SBE 50's response to each command is not shown in the examples. Review the operation of the basic sampling modes and the commands described in *Command Descriptions* before setting up your system.

Polled Sampling

Note:

Sea-Bird ships the SBE 50 with AutoRun=N (it will not automatically start sampling when power is applied). If you previously sent AutoRun=Y, send Stop to stop sampling. You will then be able to send commands to change the SBE 50 setup. On command (**TS**), the SBE 50 takes one sample of data and transmits the data real-time.

Example: Polled Sampling (user input in bold).

Set up to output pressure in decibars, with two digits after the decimal point. Set up to wait for a command upon power up. Verify setup with status (**DS**) command. Remove power.

(Apply power, then click Connect on Toolbar.) S>OUTPUTFORMAT=2 S>DECIMALS=2 S>AUTORUN=N S>DS (to verify setup) (Remove power.)

When ready to take a sample: (Apply power, then click Connect on Toolbar.) S>**TS**

Autonomous Sampling

The SBE 50 samples data at 16 Hz. The output rate is determined by NAvg:

Output rate = 16 samples per second / NAvg

AutoRun defines how sampling starts:

Note:

Sea-Bird ships the SBE 50 with AutoRun=N (it will not automatically start sampling when power is applied). If you send AutoRun=Y, turn power off and then on again to start sampling, or send Start.

AutoRun=	To Start Autonomous Sampling	To Stop Autonomous Sampling
Y	Turn on power.	 Turn off power, or Send Stop or type Ctrl Z. (Stopping this
N	Send Start.	way allows you to then send commands to check/modify the setup.)

Note:

When using the SBE 50 with an AFM and an SBE 32 Carousel Water Sampler, or using the SBE 50 with an SBE 55 ECO Water Sampler, the following setup is required:

- AutoRun=Y
- OutputFormat=7
- NAvg=16

SeatermAF automatically sets these parameters in the SBE 50 to the appropriate values when you arm the AFM or ECO, and then sends **Start** to begin sampling. *Example 1: Autonomous Sampling Setup for* **AutoRun=N** (user input in bold).

Set up to output a sample every 0.5 seconds (averaging 8 samples per output sample) and output pressure in decibars, with 2 digits after the decimal point. Verify setup with status (**DS**) command. Remove power.

(Apply power, then click Connect on Toolbar.)

S>NAVG=8

S>OUTPUTFORMAT=2 S>DECIMALS=2 S>AUTORUN=N S>DS (to verify setup) (Remove power.)

When ready to begin sampling: (Apply power, then click Connect on Toolbar.) S>START

When ready to stop sampling: (Press Enter key to get S> prompt.) S>**STOP**

Example 2: Autonomous Sampling Setup for **AutoRun=Y** (user input in bold).

Set up to output a sample every 0.5 seconds (averaging 8 samples per output sample) and output pressure in decibars, with 2 digits after the decimal point. Verify setup with status (**DS**) command. Remove power.

(Apply power, then click Connect on Toolbar.)

- S>NAVG=8
- S>OUTPUTFORMAT=2 S>DECIMALS=2 S>AUTORUN=Y S>DS (to verify setup) (Remove power.)

When ready to begin sampling: (Apply power – SBE 50 automatically begins sampling.)

When ready to stop sampling: (Press Enter key to get S> prompt.) S>**STOP**

Baud Rate, Cable Length, Power, and Data Transmission Rate

On the ship, cables longer than 3 meters should be installed inside an earthed metal conduit by a qualified electrician. This minimizes the potential for external signals to disrupt communication and ensures that high voltage lines (such as the sea cable) are sufficiently protected. Cables shorter than 3 meters can be used without shielding when installing or bench testing the instrument.

Baud Rate, Cable Length, and Data Transmission Rate

The rate that real-time data can be transmitted from the SBE 50 is dependent on the amount of data to be transmitted per scan and the serial data baud rate:

Time required to transmit data =

(number of characters * 10 bits/character) / baud rate

where

Notes:

NAvg=.

Baud rate is set with Baud=.Real-time output rate is set with

See Command Descriptions in this

• Output format is set with OutputFormat=.

section for command details.

Number of characters is dependent on the included data and output format (see *Data Output Formats* in this section). Add 2 to the number of characters shown in the output format, to account for the carriage return and line feed at the end of each scan. For decimal output (**OutputFormat=1** through **6**), include decimal points, commas, and spaces when counting the number of characters.

Time required to transmit data must be less than the real-time output rate.

The length of cable that the SBE 50 can drive to transmit real-time data is also dependent on the baud rate. The allowable combinations are:

Maximum Cable Length (meters)	Maximum Baud Rate
1600	600
800	1200
400	2400
200	4800
100	9600
50	19200
25	38400

Example - What is the fastest rate you can transmit real-time data over 800 m with **OutputFormat=0** (raw A/D pressure counts and pressure temperature voltage)?

With 800 meters of cable, the SBE 50 requires a baud rate of 1200 or less.

Number of characters for OutputFormat=0 (from Data Output Formats) =

6 (P) + 2 (comma & space) + 6 (P temperature voltage) + 2 (carriage return & line feed) = 16

Time required to transmit data = (16 characters * 10 bits/character) / 1200 =

0.133 seconds > 0.0625 seconds (16 Hz, maximum sampling rate).

Therefore, set NAvg=3, averaging 3 measurements/sample and transmitting 1 sample every 0.1875 seconds.

Note: Commor	wire resistances:
Gauge	Resistance (ohms/foot)
12	0.0016
14	0.0025
16	0.0040
18	0.0064
19	0.0081
20	0.0107
22	0.0162
24	0.0257
26	0.0410
28	0.0653

Power and Cable Length

There are two issues to consider:

- Limiting the IR loss to 1 volt; higher IR loss will prevent the instrument from transmitting real-time data because of the difference in ground potential.
- Supplying enough power at the power source so that sufficient power is available at the instrument after considering IR loss.

Looking at each issue separately:

Limiting IR Loss to 1 Volt

The limit to cable length is typically reached when the maximum current times the power common wire resistance is more than 1 volt, because the difference in ground potential of the SBE 50 and ground controller prevents the SBE 50 from transmitting real-time data.

 $V_{limit} = 1 \text{ volt} = IR_{limit}$ Where I is the current required by the SBE 50. Maximum cable length = R_{limit} / wire resistance per foot

Example 1 - For 8 volt power and 18 gauge wire, what is maximum distance to transmit power to SBE 50? Is this controlling factor for maximum cable length if wanting to transmit at 1200 baud?

With 8 volt power, the SBE 50 draws 51 milliamps (from Specifications in Section 2: Description of SBE 50). $R_{limit} = V_{limit} / I = 1 \text{ volt} / 0.051 \text{ amps} = 19.6 \text{ ohms}$

For 18 gauge wire, resistance is 0.0064 ohms/foot.

Therefore, maximum cable length is 19.6 ohms / 0.0064 ohms/foot = 3063 feet = 934 meters.

Note that 934 meters > 800 meters (maximum distance SBE 50 can transmit data at 1200 baud), so IR loss is not controlling factor for this example.

Example 2 - Same as above, but there are 4 SBE 50s powered from the same power supply.

With 8 volt power, the SBE 50 draws 51 milliamps (from Specifications in Section 2: Description of SBE 50). R limit = V limit / I = 1 volt / (0.051 amps * 4 SBE 50s) = 4.9 ohms For 18 gauge wire, resistance is 0.0064 ohms/foot.

Therefore, maximum cable length is 4.9 ohms / 0.0064 ohms/foot = 764 feet = 233 meters (this is cable length to SBE 50 furthest from power source).

Supplying Enough Power to SBE 50

Another consideration in determining maximum cable length is supplying enough power at the power source so that sufficient voltage is available, after the IR loss in the cable, to power the SBE 50.

Example 1 - For 18 gauge wire, what is maximum distance to transmit power to SBE 50 if using 8.5 volt power source to supply power? Is this controlling factor for maximum cable length if wanting to transmit at 1200 baud?

SBE 50's input power specification is 8 - 30 volts. Therefore, a 0.5 volt IR drop (8.5 volts - 8 volts) would still provide enough power to SBE 50.

0.5 volts = (.051 amps) * (0.0064 ohms/foot * cable length)V = IRCable length = 1530 ft = 466 metersNote that 466 meters < 800 meters (maximum distance SBE 50 can transmit data at 1200 baud), so IR drop in power is controlling factor for this example. Using a higher voltage power supply or a different wire gauge (12 or 14 gauge) would provide sufficient power at SBE 50 to allow 800 meter cable length.

Example 2 - Same as above, but there are 4 SBE 50s powered from the same power supply.

V = IR 0.5 volts = (.051 amps * 4 SBE 50s) *(0.0064 ohms/foot * cable length) Cable length = 382 ft = 116 meters (this is cable length to SBE 50 furthest from power source).

Optional Voltage Output Interface

The optional voltage output interface provides voltage output proportional to pressure or depth, in addition to the RS-232C output. If the voltage output feature is installed (6-pin connector in place of standard 4-pin connector), the SBE 50 can output 0.05 - 4.95 volts.

- See Dimensions and End Cap Connector in Section 2: Description of SBE 50 for voltage pinouts.
- See Optional 0-5 VDC Voltage Output Commands in Command Descriptions to enable / disable voltage output (OutputVolt=), and to set up the voltage scale factor (VRange=) to meet your deployment needs.

Command Descriptions

This section describes commands and provides sample outputs. See Appendix III: Command Summary for a summarized command list.

When entering commands:

- Input commands to the SBE 50 in upper or lower case letters and register commands by pressing the Enter key.
- The SBE 50 sends ? CMD if an invalid command is entered.
- If the system does not return an S> prompt after executing a command, press the Enter key to get the S> prompt.
- If the SBE 50 is transmitting data and you want to stop it, send Stop or type Ctrl Z. Press the Enter key to get the S> prompt. Alternatively, remove power.
- The SBE 50 responds only to Stop while sampling autonomously.

Entries made with the commands are permanently stored in the SBE 50 and remain in effect until you change them.

Status Command	
DS	Display operating status and setup parameters. Equivalent to Status on Toolbar. List below includes, where applicable, command used to modify parameter.
	 firmware version and serial number number of scans to average [NAvg=] start autonomous sampling automatically when power applied [AutoRun=]? Output format [OutputFormat=], digits after decimal point [Decimals=] Optional 0-5 VDC voltage output [OutputVolt=]? Note: This line only appears if enabled.
Example: (user input in bold)	
<pre>S>DS SBE50 V 1.0b SERIAL NO. 0011 number of scans to average = 1 start sampling on power up = no output format = depth, salt, meters, 2 decimals</pre>	[NAvg=] [AutoRun=] [OutputFormat=, Decimals=]
voltage output enabled	[OutputFormat_, Decimais_] [OutputVolt=]

	Setup Commands	
Note: The SBE 50's baud rate (set with Baud=) must be the same as	Baud=x	x = baud rate (600, 1200, 2400, 4800, 9600, 19200, or 38400). Default 9600.
SEATERM's baud rate (set in the Configure menu).	OutputFormat=x	Data output format (applies to both autonomous and polled sampling):
Nataa		x=0 : Raw pressure A/D counts and pressure temperature voltage.
 Notes: See Data Output Formats for details. If you change the data output format: Re-enter POffset= with a value that corresponds to the selected units. See Calibration Coefficients Commands below. If the optional voltage output is 		x=1 : Pressure in psia. This is required format for SBE 50 interfacing with SBE 16 <i>plus</i> , 16 <i>plus</i> V2, or 16 <i>plus</i> -IM V2, and using Sea-Bird software for real-time data acquisition and/or processing.
		x=2 (default): Pressure in decibars.
enabled, re-enter VRange= with a value that corresponds to the		x=3 : Depth in meters for salt water.
selected units. See Optional 0-5 VDC Voltage Output		x=4 : Depth in feet for salt water.
<i>Commands</i> below.When using the SBE 50 with an		x=5 : Depth in meters for fresh water.
AFM for control of an SBE 32		x=6 : Depth in feet for fresh water.
Carousel Water Sampler, or using the SBE 50 with an SBE 55 ECO Water Sampler, SeatermAF automatically sets OutputFormat=7 , NAvg=16 , and AutoRun=Y in the SBE 50 when you arm the AFM or ECO. SeatermAF then automatically sends Start to begin sampling.		x=7 : Pressure (decibars) and scan number in hexadecimal form. OutputFormat=7 required for interfacing SBE 50 with AFN and SBE 32 Carousel Water Sampler, or interfacing SBE 50 with SBE 55 ECO Water Sampler.
	Decimals=x	x= number of digits to output to right of decimal point. Applicable forOutputFormat=1 through 6.
	Examples: C	Dutput formats and digits
	• Output 553438	Format=0 , 1.5971

- OutputFormat=1 and Decimals=2 291.62
- OutputFormat=3 and Decimals=1 189.3

Note:

The input value for **Latitude=** is listed with the calibration coefficients in response to **DCal**.

Latitude=x

x= latitude (degrees) to use in salt water depth calculation. Applicable forOutputFormat=3 and 4.

- - - -

Optional 0-5 VDC Voltage Output Commands

If the optional voltage output feature is installed (6-pin connector in place of standard 4-pin connector), the SBE 50 can output 0.05 - 4.95 volts in addition to the RS-232 data. Use **OutputVolt=** to enable / disable output voltage.

OutputVolt=x	x=Y : Optional voltage output enabled.	
	x=N : Optional voltage output not enabled.	
VRange=x	x = Voltage output scale factor.	
	Output voltage = 0.1 + (4.8 * P / VRange) where	
	P = computed pressure or depth in units set with OutputFormat= .	
	If the computed voltage is less than 0.05, it is set to 0.05. If it is greater than 4.95, it is set to 4.95.	

Examples: SBE 50 with 100 meter pressure sensor.

Start by rearranging equation to solve for Vrange: **VRange** = (4.8 * P) / (Output Voltage - 0.1)

Example 1: **OutputFormat=2** (decibars), and you want to set up SBE 50 so that output is at its maximum (4.95 volts) when pressure is 100 decibars.

Substituting P = 100 db and Output voltage = 4.95 volts into equation, **VRange** = (4.8 * 100) / (4.95 - 0.1) = **98.97**. Enter this value into SBE 50.

Example 2: **OutputFormat=4** (feet, in salt water), and you want to set up SBE 50 so that output is at its maximum (4.95 volts) when depth is 328 feet.

Substituting P = 328 feet and Output voltage = 4.95 volts into equation, **VRange** = (4.8 * 328) / (4.95 - 0.1) = 324.62. Enter this value into SBE 50.

 Notes: The scale factor calculated using VCal is independent of the scale factor entered with VRange=. The values for VRange=, VOffset=, and VScale= are listed with the calibration coefficients in response to DCal. 	VCal	 Test used to internally calculate SBE 50's D/A converter offset and scale factor (VOffset= and VScale= respectively). This test requires use of a voltmeter. When VCal is sent: 1. SBE 50 outputs 0.2 volts, and then prompts for voltmeter reading. Enter voltmeter reading (in volts) and press Enter. 2. SBE 50 outputs 4.8 volts, and then prompts for voltmeter reading. Enter voltmeter reading. Enter voltmeter reading (in volts) and press Enter. 2. SBE 50 outputs 4.8 volts, and then prompts for voltmeter reading. Enter voltmeter reading. Enter voltmeter reading (in volts) and press Enter. Based on the entered voltmeter readings, the SBE 50 calculates the D/A offset and scale factor, and stores them in EEPROM.
	Volt=x	Set output voltage to x (volts). Use this command, and a voltmeter, to verify that the D/A converter offset and scale factor calculated with VCal are correct. If VCal was performed properly, the voltmeter reading should equal x .

Note:

If you change the data output format (**OutputFormat=**), re-enter **VRange=** with a value that corresponds to the selected units.

Sampling Commands

Sampling commands configure the SBE 50's response upon waking up, and direct the SBE 50 to sample data once or continuously.

Note: After you send AutoRun=Y , to start autonomous sampling	AutoRun=x	x=N : Wait for a command when power is applied. Default.
immediately: • turn power off and then on again, or • send Start		 x=Y: Start autonomous sampling at 16 Hz automatically when power is applied. SBE 50 will continue sampling until power is removed or Stop is sent. AutoRun=Y required for interfacing
Note: When using the SBE 50 with an AFM for control of an SBE 32 Carousel Water Sampler, or using the SBE 50		SBE 50 with AFM and SBE 32 Carousel Water Sampler, or interfacing SBE 50 with SBE 55 ECO Water Sampler.
with an SBE 55 ECO Water Sampler, SeatermAF automatically sets AutoRun=Y, NAvg=16, and OutputFormat=7 in the SBE 50 when you arm the AFM or ECO. SeatermAF then automatically sends Start to begin sampling.	NAvg=x	x = number of samples to average for autonomous sampling (range 1 to 255; default 1). SBE 50 samples at 16 Hz (every 0.0625 seconds), averages NAvg samples, and transmits the averaged data. NAvg=16 required for interfacing SBE 50 with AFM and SBE 32 Carousel Water Sampler, or interfacing SBE 50 with SBE 55 ECO Water Sampler.
	Start	 Start autonomous sampling at 16 Hz. Applicable if: AutoRun=N, or AutoRun=Y and you previously sent Stop to stop sampling.
 Notes: You may need to send Stop several times to get the SBE 50 to respond. You can also stop autonomous sampling by typing Ctrl Z or removing power. 	Stop	Stop autonomous sampling. Press Enter key to get an S> prompt before entering this command. Once sampling is stopped and SEATERM shows S> prompt, you can send desired commands (setup, status, etc.) to SBE 50.
	TS	Take one sample and transmit data.

DCal

Calibration Coefficients Commands

Notes:

- Date shown is when calibration was performed. Calibration coefficients are initially factory-set and should agree with Calibration Certificate shipped with SBE 50.
- See individual Coefficient Commands below for definitions of the data in the example.

Display calibration coefficients. Equivalent to Coefficients on Toolbar.

Example: (user input in bold).
S>DCAL
SBE 50 V 1.0b SERIAL NO. 0011
pressure S/N 2357, range = 3000 psia: 14-jul-01
PA0 = -9.800510e+00
PA1 = 8.881598e-03
PA2 = -1.574956e-10
PTCA0 = 5.204502e+05
PTCA1 = -7.501161e+01
PTCA2 = 1.016521e+00
PTCB0 = 2.483463e+01
PTCB1 = -7.500000e-05
PTCB2 = 0.000000e+00
PTEMPA0 = -6.530763e+01
PTEMPA1 = 5.171258e+01
PTEMPA2 = -2.201600e-01
POFFSET = 0.000000e+00
VOFFSET = -8.717179e-06
VSCALE = 1.001546e+00
VRANGE = 2.000000e+03
Latitude = 45.0

The individual Coefficient Commands listed below are used to modify a particular coefficient or date:

Note: F = floating point number S = string with no spaces	PSN=S PRange=F PCalDate=S PA0=F	S=Pressure sensor serial number F=Full scale pressure range in psia S=Pressure calibration date F=Pressure A0
	PA1=F PA2=F	F=Pressure A1 F=Pressure A2
	PTCA0=F	F=Pressure ptca0
	PTCA1=F PTCA2=F	F=Pressure ptca1 F=Pressure ptca2
	PTCB0=F PTCB1=F	F=Pressure ptcb0
	PTCB2=F	F=Pressure ptcb1 F=Pressure ptcb2
	PTempA0=F PTempA1=F	F=Pressure ptempa0 F=Pressure ptempa1
Г	PTempA2=F	F=Pressure ptempa2
Note: If you change the data output format (OutputFormat=), re-enter POffset=	POffset=F	F=Pressure offset; added to computed pressure in units determined with OutputFormat =
with a value that corresponds to the selected units.	VOffset=, VScale=, VRange=	See Optional 0-5 VDC Voltage Output Commands above.
	Latitude=	See Setup Commands above.

Data Output Formats

Notes:

- Each scan of data ends with a carriage return and line feed.
- OutputFormat=0 through 6: data is in decimal form.
- OutputFormat=1 through 6: Leading zeros are suppressed, except for one zero to the left of the decimal point. Decimals= defines the number of digits to the right of the decimal point.
- **OutputFormat=1** (psia) is required format for an SBE 50 interfacing with an SBE 16*plus*, 16*plus* V2, or 16*plus*-IM V2 and using Sea-Bird software for real-time data acquisition and/or processing.

Note:

When using the SBE 50 with an AFM for control of an SBE 32 Carousel Water Sampler, or using the SBE 50 with an SBE 55 ECO Water Sampler, SeatermAF automatically sets **OutputFormat=7**, **NAvg=16**, and **AutoRun=Y** in the SBE 50 when you arm the AFM or ECO. SeatermAF then automatically sends **Start** to begin sampling. • **OutputFormat=0** - Raw A/D pressure counts and pressure temperature voltage:

pppppp, v.vvvv

where pppppp = pressure counts v.vvvv = pressure temperature voltage The pressure sensor is an absolute sensor, so its raw output includes the effect of atmospheric pressure (14.7 psi). See the SBE 50's Calibration Sheet for the equation to convert the pressure counts and pressure temperature voltage to pressure.

• **OutputFormat=1** (psia) and **2** (decibars) - Pressure: When **OutputFormat=1**, the SBE 50 outputs pressure in psia, which includes the effect of atmospheric pressure (i.e., output pressure at the surface is 14.7 psia).

When **OutputFormat=2**, the SBE 50 outputs pressure in decibars, relative to the ocean surface (i.e., output pressure at the surface is 0 decibars).

As shown on the Calibration Sheet, Sea-Bird's calibration (and resulting calibration coefficients) is in terms of psia. The SBE 50 uses the following equation to convert psia to decibars:

pressure (db) = [pressure (psia) - 14.7] * 0.689476; Pressure in decibars is the value required by the UNESCO formulas for computation of salinity, density, and other derived variables.

- **OutputFormat=3** (meters) and **4** (feet) Depth for salt water: For oceanic applications, it is assumed that salinity is 35 PSU and temperature is 0 °C, and the water compressibility (with accompanying density variation) is taken into account in the depth calculation. This is the method recommended in UNESCO Technical Paper No. 44 and is a logical approach, because most of the deep-ocean water column approximates these salinity and temperature values. Since pressure is also proportional to gravity and the major variability in gravity depends on latitude, the input latitude (**Latitude=**) is used to estimate the magnitude of the local gravity.
- **OutputFormat=5** (meters) and **6** (feet) Depth for fresh water: For fresh water applications, compressibility is not significant in the shallow depths encountered and is ignored, as is the latitude-dependent gravity variation. Fresh water density is assumed to be 1 gm/cm for the calculation of depth.
- **OutputFormat=7** Pressure and scan number in hexadecimal: ppppsssss

where pressure [decibars] = pppp - 100 ssssss = scan number
OutputFormat=7 required for communicating with the Auto Fire
Module (AFM) and SBE 32 Carousel Water Sampler or the SBE 55
ECO Water Sampler. This allows the AFM or ECO to use the SBE 50's pressure data to determine when to close water bottles.

Example: SBE 50 used with AFM and Carousel, example scan = ppppsssss = 00C80001F0

- Pressure = pppp = 00C8 (200 decimal); pressure (decibars) = 200 - 100 = 100 decibars
- Scan number = ssssss = 0001F0 (496 decimal); scan number = 496


SBE 50

Deployment

CAUTION:

Do not use WD-40 or other petroleumbased lubricants, as they will damage the connector.





- 1. Install the I/O cable on the SBE 50 end cap:
 - A. Lightly lubricate the inside of the cable connector with silicone grease (DC-4 or equivalent).
 - B. **Standard Connector** Install the cable connector, aligning the raised bump on the side of the connector with the large pin (pin 1 ground) on the SBE 50. Remove any trapped air by *burping* or gently squeezing the connector near the top and moving your fingers toward the end cap. **OR**

MCBH Connector – Install the cable connector, aligning the pins.

- C. Place the locking sleeve over the connector. Tighten the locking sleeve finger tight only. Do not overtighten the locking sleeve and do not use a wrench or pliers.
- 2. Connect the I/O cable to the computer/controller.
- 3. Connect the I/O cable's red (+) and black (-) wires to a power supply (8-30 VDC).
- 4. Verify that the hardware and external fittings are secure.
- 5. Immediately prior to deployment, apply power.
 - If AutoRun=Y, the SBE 50 will begin sampling.
 - If AutoRun=N, send Start to begin autonomous (continuous) sampling or send TS to take a single sample.
- 6. If the SBE 50 is connected directly to the computer and you will be using SEATERM to view real-time data, press Capture on the Toolbar to save the data to a .cap file.

The SBE 50 is ready to go into the water.

Recovery

WARNING!

If the SBE 50 stops working while underwater, is unresponsive to commands, or shows other signs of flooding or damage, carefully secure it away from people until you have determined that abnormal internal pressure does not exist or has been relieved. Pressure housings may flood under pressure due to dirty or damaged o-rings, or other failed seals. When a sealed pressure housing floods at great depths and is subsequently raised to the surface, water may be trapped at the pressure at which it entered the housing, presenting a danger if the housing is opened before relieving the internal pressure. Instances of such flooding are rare. However, a housing that floods at 5000 meters depth holds an internal pressure of more than 7000 psia, and has the potential to eject the end cap with lethal force. A housing that floods at 50 meters holds an internal pressure of more then 85 psia; this force could still cause injury. If you suspect the SBE 50 is flooded,

point it in a safe direction away from people, and loosen the bulkhead connector very slowly, at least 1 turn. This opens an o-ring seal under the connector. Look for signs of internal pressure (hissing or water leak). If internal pressure is detected, let it bleed off slowly past the connector o-ring. Then, you can safely remove the end cap. Rinse the outside of the SBE 50 with fresh water.

Section 5: Routine Maintenance and Calibration

This section reviews corrosion precautions, connector mating and maintenance, pressure sensor maintenance, and sensor calibration. The accuracy of the SBE 50 is sustained by establishing proper handling practices.

Corrosion Precautions

Rinse the SBE 50 with fresh water after use and prior to storage.

All exposed metal is titanium; other materials are plastic. No corrosion precautions are required, but direct electrical connection of titanium to dissimilar metal hardware should be avoided.

Connector Mating and Maintenance

Note:

See Application Note 57: Connector Care and Cable Installation.

CAUTION:

Do not use WD-40 or other petroleum-based lubricants, as they will damage the connector.

Clean and inspect the connectors, cable, and dummy plug before every deployment and as part of your yearly equipment maintenance. Inspect connectors that are unmated for signs of corrosion product around the pins, and for cuts, nicks or other flaws that may compromise the seal.

When remating:

- 1. Lightly lubricate the inside of the cable connector with silicone grease (DC-4 or equivalent).
- Standard Connector Install the cable connector, aligning the raised bump on the side of the cable connector with the large pin (pin 1 ground) on the SBE 50. Remove any trapped air by *burping* or gently squeezing the connector near the top and moving your fingers toward the end cap. OR
 MCBH Connector Install the cable connector, aligning the pins.
- 3. Place the locking sleeve over the cable connector. Tighten the locking sleeve finger tight only. Do not overtighten the locking sleeve and do not use a wrench or pliers.

Verify that a cable is installed before deployment.

Pressure Sensor Maintenance



CAUTION: Do not put a brush or any object in the pressure port. Doing so may damage or break the pressure sensor. The pressure port plug has a small vent hole to allow hydrostatic pressure to be transmitted to the pressure sensor inside the instrument, while providing protection for the pressure sensor, keeping most particles and debris out of the pressure port.

Periodically (approximately once a year) inspect the pressure port to remove any particles, debris, etc:

- 1. Unscrew the pressure port plug from the pressure port.
- 2. Rinse the pressure port with warm, de-ionized water to remove any particles, debris, etc.
- 3. Replace the pressure port plug.

Sensor Calibration



Note:

The SBE 50's pressure sensor is an absolute sensor, so its raw output includes the effect of atmospheric pressure (14.7 psi). As shown on the Calibration Sheet, Sea-Bird's calibration (and resulting calibration coefficients) is in terms of psia. However, when outputting pressure in engineering units, the SBE 50 outputs pressure relative to the ocean surface (i.e., at the surface the output pressure is 0 decibars). The SBE 50 uses the following equation to convert psia to decibars: Pressure (db) =[pressure (psia) - 14.7] * 0.689476

Sea-Bird sensors are calibrated by subjecting them to known physical conditions and measuring the sensor responses. Coefficients are then computed, which may be used with appropriate algorithms to obtain engineering units. The SBE 50's pressure sensor is supplied fully calibrated, with coefficients printed on the Calibration Certificate. These coefficients have been stored in the SBE 50's EEPROM.

The strain-gauge pressure sensor is a mechanical diaphragm type, with an initial static error band of 0.05%. Consequently, the sensor is capable of meeting SBE 50's 0.10% error specification with some allowance for aging and ambient-temperature induced drift.

Pressure sensors show most of their error as a linear offset from zero. A technique is provided below for making small corrections to the pressure sensor calibration using the *offset* (**POffset=**) calibration coefficient term by comparing SBE 50 pressure output to readings from a barometer.

Allow the SBE 50 to equilibrate (with power on) in a reasonably constant temperature environment for at least 5 hours before starting. Pressure sensors exhibit a transient change in their output in response to changes in their environmental temperature. Sea-Bird instruments are constructed to minimize this by thermally decoupling the sensor from the body of the instrument. However, there is still some residual effect; allowing the SBE 50 to equilibrate before starting will provide the most accurate calibration correction.

- 1. Place the SBE 50 in the orientation it will have when deployed.
- 2. In SEATERM:
 - A. Set the pressure offset to 0.0 (**POffset=0**).
 - B. Send TS a number of times to measure and transmit SBE 50 pressure.
- Compare the SBE 50 output to the reading from a good barometer at the same elevation as the SBE 50's pressure sensor. Convert the barometer reading to units consistent with **OutputFormat=**. Calculate *offset* = barometer reading SBE 50 reading
- 4. Enter the offset (positive or negative) in the SBE 50's EEPROM, using **POffset=** in SEATERM. **The offset must be entered in units consistent** with **OutputFormat=** (for example, if the output format is decibars, enter an offset in decibars).

Offset Correction Example

Absolute pressure measured by a barometer is 1010.50 mbar. Pressure displayed from SBE 50 with output in decibars (**OutputFormat=2**) is -2.5 dbars. Convert barometer reading to dbars using the relationship: mbar * 0.01 = dbar Barometer reading = 1010.50 mbar * 0.01 = 10.1050 dbar The SBE 50's internal calculations and our processing software output gage pressure, using an assumed value of 14.7 psi for atmospheric pressure. Convert SBE 50 reading from gage to absolute by adding 14.7 psia to the SBE 50's output: -2.5 dbars + (14.7 psi * 0.689476 dbar/psia) = -2.5 + 10.13 = 7.635 dbars Offset = 10.1050 - 7.635 = + 2.47 dbars Enter offset in SBE 50 with **POffset=**.

For demanding applications, or where the sensor's air ambient pressure response has changed significantly, calibration using a dead-weight generator is recommended. The pressure sensor port uses a 7/16-20 straight thread for mechanical connection to the pressure source. Use a fitting that has an O-ring tapered seal, such as Swagelok-200-1-4ST, which conforms to MS16142 boss.

Section 6: Troubleshooting

This section reviews common problems in operating the SBE 50, and provides the most likely causes and solutions.

Problem 1: Unable to Communicate with SBE 50

The S> prompt indicates that communications between the SBE 50 and computer have been established. Before proceeding with troubleshooting, attempt to establish communications again by clicking Connect on SEATERM's toolbar or hitting the Enter key several times.

Cause/Solution 1: The I/O cable connection may be loose. Check the cabling between the SBE 50 and computer for a loose connection.

Cause/Solution 2: The instrument type and/or its communication settings may not have been entered correctly in SEATERM. Select *SBE 50* in the Configure menu and verify the settings in the Configuration Options dialog box. The settings should match those on the instrument Configuration Sheet.

Cause/Solution 3: The I/O cable may not be the correct one. The I/O cable supplied with the SBE 50 permits connection to DB-9P input connectors, i.e., standard RS-232 interfaces.

- SBE 50 Pin 1 goes to DB-9 pin 5 (ground)
- SBE 50 pin 2 goes to DB-9 pin 3
- SBE 50 pin 3 goes to DB-9 pin 2

See Section 2: Description of SBE 50 for the connector details.

Problem 2: Unreasonable Data

The symptom of this problem is data that contains unreasonable values (for example, values that are outside the expected range of the data).

Cause/Solution 1: Data with unreasonable values may be caused by incorrect calibration coefficients in the SBE 50's EEPROM (if looking at data output in converted format, **OutputFormat = 1** through **7**). Verify the calibration coefficients in EEPROM, using **DCal**.

Glossary

SBE 50 – High-accuracy Digital Oceanographic Pressure Sensor.

PCB – Printed Circuit Board.

Scan – One data sample.

SEASOFT V2– Sea-Bird's complete Win 2000/XP software package, which includes software for communication, real-time data acquisition, and data analysis and display. SEASOFT V2 includes **SEATERM**, **SeatermAF**, Seasave, SBE Data Processing, and Plot39. Only SEATERM and SeatermAF are applicable to the SBE 50.

SEATERM – Sea-Bird's Win 95/98/NT/2000/XP terminal program software used to communicate with the SBE 50.

SeatermAF – Sea-Bird's Win 95/98/NT/2000/XP terminal program software used to communicate with the SBE 50 when it is used with an SBE 32 Carousel Water Sampler and Auto Fire Module (AFM), or with an SBE 55 ECO Water Sampler. See the AFM manual or the SBE 55 ECO manual for details.

Appendix I: Functional Description

Sensor

Note:

Pressure ranges are expressed in meters of deployment depth capability.

The SBE 50's pressure sensor, developed by Druck, Inc., has a superior design that is entirely different from conventional 'silicon' types in which the deflection of a metallic diaphragm is detected by epoxy-bonded silicon strain gauges. The Druck sensor employs a micro-machined *silicon diaphragm* into which the strain elements are implanted using semiconductor fabrication techniques. Unlike metal diaphragms, silicon's crystal structure is perfectly elastic, so the sensor is essentially free of pressure hysteresis. Compensation of the temperature influence on pressure offset and scale is performed by the SBE 50's CPU. The pressure sensor is available in the following pressure ranges: 20, 100, 350, 600, 1000, 2000, 3500, and 7000 meters.

Sensor Interface

Strain-gauge pressure is acquired by applying an AC excitation to the pressure bridge. A 24-bit A/D converter digitizes the output of the bridge. AC excitation and ratiometric comparison avoids errors caused by parasitic thermocouples, offset voltages, leakage currents, and reference errors. A silicon diode embedded in the pressure bridge is used to measure the temperature of the pressure bridge. This temperature is used to perform offset and span corrections on the measured pressure signal.

Appendix II: Electronics Disassembly/Reassembly

Disassembly



└ Unscrew end cap

- 1. Wipe the outside of the end cap and housing dry, being careful to remove any water at the seam between them.
- 2. Using a wrench, unscrew the end cap.
- 3. Carefully pull the end cap and attached electronics out of the housing. Note that the electronics are electrically connected to the I/O connector.
- 4. Remove any water from the O-ring and mating surfaces inside the housing with a lint-free cloth or tissue.
- 5. Be careful to protect the O-ring from damage or contamination.



Reassembly

Note:

Before delivery, a desiccant package is placed in the housing, and the housing is filled with dry Argon gas. These measures help prevent condensation.

To ensure proper functioning:

- Install a new desiccant bag each time you open the housing. If a new bag is not available, see Application Note 71: Desiccant Use and Regeneration (drying).
- If possible, dry gas backfill each time you open the housing. If you cannot, wait at least 24 hours before redeploying, to allow the desiccant to remove any moisture from the housing.
- 1. Remove any water from the end cap's O-ring and mating surfaces in the housing with a lint-free cloth or tissue. Inspect the O-ring and mating surfaces for dirt, nicks, and cuts. Clean as necessary. Apply a light coat of O-ring lubricant (Parker Super O Lube) to the O-rings and mating surfaces.
- 2. Carefully fit the electronics into the housing until you feel resistance.
- 3. Using a wrench, screw the end cap onto the housing.

Appendix III: Command Summary

Category	Command	Description		
Status	DS	Display status and setup parameters.		
	Baud=x	x = baud rate (600, 1200, 2400, 4800, 9600, 19200, or 38400). Default 9600.		
Setup	OutputFormat= x	 x=0: raw pressure and pressure temperature. x=1: pressure in psia. x=3: depth in meters, salt water. x=5: depth in meters, fresh water. x=6: depth in feet, fresh water. x=7: pressure (decibars) and scan number in hexadecimal. 		
	Decimals=x	x = number of digits to output to right of decimal point (for OutputFormat=1 through 6).		
	Latitude=x	$\mathbf{x} = $ latitude (deg) for depth calculation (for OutputFormat=3 and 4).		
	OutputVolt=x	x=Y : Optional voltage output enabled. x=N : Not enabled.		
	VRange=x	x = Voltage output scale factor. Output voltage = $0.1 + (4.8 * P / VRange)$ where P = computed pressure or depth in units set with OutputFormat = If computed voltage < 0.05, it is set to 0.05. If > 4.95, it is set to 4.95.		
Optional 0-5 VDC Voltage Output	VCal	 Test used to internally calculate SBE 50's D/A converter offset and scale factor. Test requires use of a voltmeter. When VCal is sent: SBE 50 outputs 0.2 volts. Enter voltmeter reading (in volts) at prompt and press Enter. SBE 50 outputs 4.8 volts. Enter voltmeter reading (in volts) at prompt and press Enter. 		
	Volt=x	Set output voltage to \mathbf{x} (volts). Use Volt =, and a voltmeter, to verify that D/A converter offset and scale factor calculated with VCal are correct. If VCal performed properly, voltmeter reading should equal \mathbf{x} .		
	AutoRun=x	x=N: Wait for command when power applied. Default.x=Y: Start autonomous sampling automatically when power applied.		
Sampling	NAvg=x	x = number of samples to average for autonomous sampling. SBE 50 samples at 16 Hz, averages NAvg samples, and transmits averaged data. Default=1.		
Sumpling	Start	Start autonomous sampling now.		
	Stop	Stop autonomous sampling. Press Enter key to get S> before entering Stop .		
1	TS	Take 1 sample and transmit data.		
Coefficients (F=floating	DCal	Display calibration coefficients. Use individual commands below to modify a particular coefficient or date.		
point number;	PSN=S	F= Pressure sensor serial number.		
S=string with no spaces)	I Kange-F	F = Full scale pressure range in psia.		
no spaces)	PCalDate=S	S= Pressure calibration date.		
Date shown	PA0=F	F= Pressure A0.		
is when	PA1=F	F= Pressure A1.		
calibration	PA2=F	F = Pressure A2.		
was	PTCA0=F	F = Pressure ptca0.		
performed. Calibration	PTCA1=F	F= Pressure ptca1.		
coefficients	PTCA2=F	F= Pressure ptca2.		
are initially	PTCB0=F	F = Pressure ptcb0.		
factory-set	PTCB1=F	F= Pressure ptcb1.		
and should	PTCB2=F	F = Pressure ptcb2.		
agree with Calibration	PTempA0=F	F = Pressure ptempa0.		
Calibration	PTempA1=F	F = Pressure ptempa1.		
shipped with	PTempA2=F	F= Pressure ptempa2.		
SBE 50.	POffset=F	F = Pressure offset.		

Appendix IV: Replacement Parts

Part Number	Part	Application Description	Quantity in SBE 50
30931	Parker 2-214N674-70 O-ring	End cap O-ring	1
801385	4-pin RMG-4FS to 9-pin DB-9S I/O cable, 2.4 m (8 ft long)	SBE 50 (standard connector) to computer and power supply	1
801386	6-pin AG-206 to 9-pin DB-9S I/O cable, 2.4 m (8 ft long)	SBE 50 (with optional Voltage Output, standard connector) to computer and power supply	1
801206	4-pin MCIL-4FS to 9-pin DB-9S I/O cable, 2.4 m (8 ft long)	SBE 50 (wet-pluggable connector) to computer and power supply	1
801387	6-pin MCIL-6FS to 9-pin DB-9S I/O cable, 2.4 m (8 ft long)	SBE 50 (with optional Voltage Output, wet-pluggable connector) to computer and power supply	1
171888	25-pin DB-25S to 9-pin DB-9P cable adapter	For use with computer with DB-25 connector	1
171730	Double Y-cable (AFM 6-pin AG-206 and 3-pin RMG-3FS to SBE 32 6-pin AG-206 and SBE 50 4-pin RMG-4FS)	SBE 50 used with SBE 32 Carousel Water Sampler and Auto Fire Module (standard connectors on SBE 50, SBE 32, and AFM)	-
171991	Double Y-cable (AFM 6-pin MCIL-6FS and 3-pin MCIL-3FS to SBE 32 6-pin MCIL-6FS and SBE 50 4-pin MCIL-4FS)	SBE 50 used with SBE 32 Carousel Water Sampler and Auto Fire Module (wet-pluggable connectors on SBE 50, SBE 32, and AFM)	-
172259	6-pin AG-206 to 4-pin RMG-4FS, 1.2 m (4 ft)	SBE 50 used with SBE 55 ECO Water Sampler (standard connectors on SBE 50 and ECO)	-
172260	6-pin MCIL-6FS to 4-pin MCIL-4FS, 1.2 m (4 ft)	SBE 50 used with SBE 55 ECO Water Sampler (wet-pluggable connectors on SBE 50 and ECO)	-
17043	Locking sleeve	Locks cable or dummy plug in place (SBE 50 with standard connector)	1
171192	Locking sleeve	Locks cable or dummy plug in place (SBE 50 with wet- pluggable connector)	1
17046.1	4-pin RMG-4FSD dummy plug with locking sleeve	Connector protection (standard connector)	1
17047.1	6-pin AG-206 dummy plug with locking sleeve	Connector protection (SBE 50 with optional Voltage Output, standard connector)	1
171398.1	4-pin MCDC-4-F dummy plug with locking sleeve	Connector protection (wet-pluggable connector) when I/O cable not connected	1
171498.1	6-pin MCDC-6-F dummy plug with locking sleeve	Connector protection when I/O cable not connected (SBE 50 with optional Voltage Output, wet-pluggable connector)	1

Index

0

0 - 5 Volt output \cdot 18, 20

A

About Sea-Bird · 4 AFM · 7, 15, 23 Auto Fire Module · 7, 15, 23

B

Baud rate \cdot 16

С

Cable length · 16, 17 Calibration · 28 Carousel Water Sampler · 7, 15, 23 Command summary · 33 Commands calibration coefficients \cdot 22 descriptions · 18 operating \cdot 21 output format · 19 sampling $\cdot 21$ setup · 19 status · 18 voltage output \cdot 20 Communication defaults · 12 Connector · 9, 26 Corrosion precautions · 26

D

Data output format \cdot 19, 23 Data transmission rate \cdot Deployment \cdot preparing for \cdot Description \cdot Dimensions \cdot

E

ECO Water Sampler \cdot 7, 15, 23 Electronics disassembly/reassembly \cdot 32 End cap \cdot 26

F

Functional description \cdot 31

G

Glossary · 30

L

Limited liability statement $\cdot \ 2$

М

Maintenance \cdot 26 Modes \cdot *See* Sampling modes

0

Optional voltage output interface \cdot 18, 20 Output format \cdot 23

Р

Parts replacement · 35 Power · 8, 17 Pressure sensor maintenance · 27

Q

Quick start · 4

R

Recovery · 25 Replacement parts · 35

S

Sampling modes · 14 autonomous · 15 polled · 14 SBE 32 Carousel Water Sampler · 7, 15, 23 SBE 55 ECO Water Sampler · 7, 15, 23 SEASOFT · 7, 10 SEATERM · 7, 10, 11 main screen · 11 toolbar buttons · 12 SeatermAF · 7 Sensor · 8 Software · 7 Software installation · 10 Specifications · 8 System description · 6

T

Terminal program installation · 10 Testing · 10 Troubleshooting · 29

U

Unpacking SBE 50 · 5

V

Voltage output interface · 18, 20

W

Water sampler \cdot 7, 15, 23 Wiring \cdot 10