

YYYY

MM

DD

HH

MM

#

OPERATOR(S) INITIALS

2023 08 21 12 00

_STATION-

066

_SML

E. Boss

Water Collection

Number of dips: 13

Volume Collected (L)

700 ml

Protocol	Quantity, Container Storage	Replicate 1	Replicate 2	Replicate 3	Comments
SML-FC	Cryotube (2ml) 1.5ml sample + 30 µl Glute 25% LN2 (-80°C)				Strongly with viral
SML-CP	Cryotube (5ml) 3ml sample + 750 µl Glycerol RT (-80°C)				
SML-320	Cryotube (2ml) Filter 3µm PC LN2 (-80°C)				~66 ml
SML-023	Cryotube (2ml) Filter 0.22µm PC LN2 (-80°C)				~60 ml

YYYY

MM

DD

HH

MM

#

OPERATOR(S) INITIALS

_STATION-

_SML

Water Collection

Number of dips:

Volume Collected (L)

Protocol	Quantity, Container Storage	Replicate 1	Replicate 2	Replicate 3	Comments
SML-FC	Cryotube (2ml) 1.5ml sample + 30 µl Glute 25% LN2 (-80°C)	SML-FC-1	SML-FC-2		
SML-CP	Cryotube (5ml) 3ml sample + 750 µl Glycerol RT (-80°C)	SML-CS-1	SML-CS-2	SML-CS-3	
SML-320	Cryotube (2ml) Filter 3µm PC LN2 (-80°C)	SML-320-1	SML-320-2	SML-320-3	
SML-023	Cryotube (2ml) Filter 0.22µm PC LN2 (-80°C)	SML-023-1	SML-023-2	SML-023-3	



COMMENTS

1	
2	
3	
4	
5	
6	
7	