



YYYY MM DD # # #
 LOG_SAMPLES_ 2023 09 25 _STATION- 080 _W-LAB-142-1
 OPERATOR(S) Julie Paulain

Depth Replicate	S320 Cryo-5mL LN2 #1	S023 Cryo-5mL LN2 #1	Filtration Volume (Litres)	Filtration Duration (minutes)	S<02 Cryo-5mL FRG +4°C	Filtration Volume (Litres)
Z00 R01 m			[] 20L [] 50L 10 L	[] 15' [] 60' 15 min.		[] 10L [] 20L 10 L
Z00 R02 m			[] 20L [] 50L 10 L	[] 15' [] 60' 15 min.		[] 10L [] 20L 10 L
Z02 R01 m	###-Z02 S320-1	###-Z02 S023	[] 20L [] 50L L	[] 15' [] 60' min.	###-Z02 S<02	[] 10L [] 20L L
Z02 R02 m	###-Z02 S320-2	###-Z02 S023-2	[] 20L [] 50L L	[] 15' [] 60' min.	###-Z02 S<02-2	[] 10L [] 20L L
Depth Replicate	P320 Cryo-5mL LN2 #1	P023 Cryo-5mL LN2 #1	Filtration Volume (Litres)	Filtration Duration (minutes)	< 0.2 µm	
Z00 m			[] 20L [] 50L 25 L	[] 15' [] 60' 60 min.	=> Collect filtrate for SS protocols onland : VV<0.2, qPCR<0.2	
Z02 m	###-Z02 P320	###-Z02 P023	[] 20L [] 50L L	[] 15' [] 60' min.		
Depth Replicate	S320-L 15mL falcon -20°C + 10 mL Nucleoprotect	S023-L 15mL falcon -20°C + 10 mL Nucleoprotect	Filtration Volume (Litres)	Filtration Duration (minutes)		
Z00 m			[] 20L [] 50L 40 L	[] 15' [] 60' 60 min.	→	
Z02 m	###-Z02 S320-L	###-Z02 S023-L	[] 20L [] 50L L	[] 15' [] 60' min.		



Depth Replicate	COMMENTS S###
Z00 R01 m	
Z00 R02 m	
Z02 R01 m	
Z02 R02 m	
Depth	COMMENTS P###
Z00 m	
Z02 m	
Depth	COMMENTS S###L
Z00 m	2 filters per tube S320 x S023
Z02 m	



LOG_SAMPLES_ YYYY MM DD # # # _STATION- # # # _DECK-BGC

OPERATOR(S) MH

2023 09 25

0 8 0

_STATION- 0 8 0 _DECK-BGC

OPERATOR(S) MH

Depth	TOC Vial-40mL FRG +4°C	TOC Vial-40mL FRG +4°C	TOC Vial-40mL FRG +4°C		DICTA Bottle-500mL RT >10°C	SAL Bottle-125mL RT >10°C
Z00 m						
Z02 m	###-Z02 TOC-1	###-Z02 TOC-2	###-Z02 TOC-3		###-Z02 DICTA	###-Z02 SAL
	+ 150 µl HCl				+ 300 µl HgCl ₂	
Depth	CDOM/FDOM Bottle-60mL FRG +4°C	DOC Vial-40mL FRG +4°C	NUT Bottle-60mL FRZ -20°C			
Z00 R01 m						
Z00 R02 m						
Z00 R03 m						
Z02 R01 m	###-Z02 DOM-1	###-Z02 DOC-1	###-Z02 NUT-1			
Z02 R02 m	##-Z02 DOM-2	###-Z02 DOC-2	###-Z02 NUT-2			
Z02 R03 m	##-Z02 DOM-3	###-Z02 DOC-3	###-Z02 NUT-3			
	+ 150 µl HCl					



Depth		COMMENTS TOC	COMMENTS DICTA	COMMENTS SAL
Z00	m			
Z02	m			
Depth Replicate		COMMENTS CDOM/FDOM	COMMENTS DOC	COMMENTS NUT
Z00	R01			
	m			
Z00	R02			
	m			
Z00	R03			
	m			
Z02	R01			
	m			
Z02	R02			
	m			
Z02	R03			
	m			



LOG_SAMPLES_ **2023 09 25** _STATION- **080** _TARDIS-SCP
 OPERATOR(S) **MH**

Depth	PPL Falcon-50mL FRZ -20°C	Filtration Volume (Litres)	Filtration Duration (minutes)	HLB Falcon-50mL FRZ -20°C	Filtration Volume (Litres)	Filtration Duration (minutes)
Z00 R01 m		[] 1L [x] 2L L	[x] 60' [] 120' min.		[] 1L [x] 2L L	[x] 60' [] 120' min.
Z00 R02 m		[] 1L [x] 2L L	[x] 60' [] 120' min.		[] 1L [x] 2L L	[x] 60' [] 120' min.
Z00 R03 m		[] 1L [x] 2L L	[x] 60' [] 120' min.		[] 1L [x] 2L L	[x] 60' [] 120' min.
Z00 R04 m		[] 1L [x] 2L L	[x] 60' [] 120' min.		[] 1L [x] 2L L	[x] 60' [] 120' min.
Z02 R01 m	###-Z02 PPL-1	[] 1L [] 2L L	[] 60' [] 120' min.	###-Z02 HLB-1	[] 1L [] 2L L	[] 60' [] 120' min.
Z02 R02 m	###-Z02 PPL-2	[] 1L [] 2L L	[] 60' [] 120' min.	###-Z02 HLB-2	[] 1L [] 2L L	[] 60' [] 120' min.
Z02 R03 m	###-Z02 PPL-3	[] 1L [] 2L L	[] 60' [] 120' min.	###-Z02 HLB-3	[] 1L [] 2L L	[] 60' [] 120' min.
Z02 R04 m	###-Z02 PPL-4	[] 1L [] 2L L	[] 60' [] 120' min.	###-Z02 HLB-4	[] 1L [] 2L L	[] 60' [] 120' min.
Depth Replicate	MB320 50mL-Falcon FRZ -20°C	Filtration Volume (Litres)	Filtration Duration (minutes)	MB033 50mL-Falcon FRZ -20°C	Filtration Volume (Litres)	Filtration Duration (minutes)
Z00 m		[x] 16L L	[] 30' min. 20		[x] 16L L	[] 30' min. 20
Z02 m	###-Z02 MB320	[] 16L L	[] 30' min.	###-Z02 MB033	[] 16L L	[] 30' min.



Depth		COMMENTS
Z00	R01 m	
Z00	R02 m	
Z00	R03 m	
Z00	R04 m	
Z02	R01 m	
Z02	R02 m	
Z02	R03 m	
Z02	R04 m	



LOG_SAMPLES_ YYYY MM DD # # # _STATION- # # # _S-LAB-OTHER

OPERATOR(S) Sc

Depth Replicates	HC Cryo-5mL LN2 #1	HC-G Cryo-5mL LN2 #1	CP-G Cryo-5mL LN2 #1	SG Cryo-5mL LN2 #1	FC-P Cryo-2mL LN2 #3	FC-G Cryo-2mL LN2 #3
Z00 R01 m						
Z00 R02 m						
Z00 R03 m				Glycine-betaine prealiquot at 4°C	PFA prealiquot at -20°C	Glutaraldehyde prealiquot at -20°C
Z00 R04 m			Glycerol prealiquot - RT			
Z00 R05 m				DGAS* 12 mL exetainer +4°C	DGAS* 12 mL exetainer +4°C	DGAS* 12 mL exetainer +4°C
Z00 R06 m						
Z00 R07 m						
Z00 R08 m						
Prealiquot	No prealiquot	Glycerol prealiquot - RT		* + 100 µL ZnCl ²	* + 100 µL ZnCl ²	* + 100 µL ZnCl ²
Depth Replicate	eDNA Watera capsule RT	Filtration Volume (Litres)	Filtration Duration (minutes)	+ 50 mL of buffer	< 0.45 µm	
Z00 m		[] 30L	33 min.	=> Collect filtrate for SS protocol onland : V<0.45		
Z02 m	###-Z02 eDNA	[] 30L	[] 30 min.			



Depth Replicate	COMMENTS
Z00 R01 m	
Z00 R02 m	
Z00 R03 m	
Z00 R04 m	
Z00 R05 m	
Z00 R06 m	
Z00 R07 m	
Z00 R08 m	



LOG_SAMPLES_ **2023** | **09** | **25** | _STATION- **0** | **8** | **0** | _S-LAB-25-1

OPERATOR(S) | **S. Z.** |

Depth	Turbidimeter (FNU)	PM control (EVERY TWO STATIONS)	Filtration Volume (mL)	N° filtres + weight (mg)		
Z00 m	1. 3.76 2. 4.29 3. 4.09	###-Z00 PM-CTRL	mL	N°: Weight:		
Z02 m	1. 2. 3.	TRIPPLICATES ONCE A MONTH FOR HP	HP Cryo-2mL LN2 #2	Filtration Volume (mL) Filtration Duration (minutes)		
Depth	PA Petridish FRZ -20°C	Filtration Volume (mL)	Z00 R01 m	[] 30' [] 40' max 15 min		
Z00 m		635 mL	Z00 R02 m	###-Z00 HP-2 mL [] 30' [] 40' max min		
Z02 m	###-Z02 PA	mL	Z00 R03 m	###-Z00 HP-3 mL [] 30' [] 40' max min		
Depth	PM Petridish FRZ -20°C	Filtration Volume (mL)	N° filtre + weight (mg)	FOI Petridish FRZ -20°C	Filtration Volume (mL)	N° filtre + weight (mg)
Z00 R01 m		635 mL	N°: XX091 Weight: 36,990		400 mL	N°: XX163 Weight: 37,669
Z00 R02 m		635 mL	N°: XX400 Weight: 37,449		400 mL	N°: XX161 Weight: 36,899
Z00 R03 m		400 mL	N°: XX316 Weight: 37,297		400 mL	N°: XX010 Weight: 36,826
Z02 R01 m	###-Z02 PM-1	mL	N°: XX Weight:	###-Z02 FOI-1	mL	N°: Weight:
Z02 R02 m	###-Z02 PM-2	mL	N°: Weight:	###-Z02 FOI-2	mL	N°: Weight:
Z02 R03 m	###-Z02 PM-3	mL	N°: Weight:	###-Z02 FOI-3	mL	N°: Weight:



Depth Replicate		COMMENTS PM	COMMENTS FOI
Z00	R01 m		
Z00	R02 m		
Z00	R03 m		
Z02	R01 m		
Z02	R02 m		
Z02	R03 m		
Depth Replicate		COMMENTS PA - HP	
Z00	m		
Z02	m		



LOG-SAMPLES_ YYYY MM DD # # # _STATION- _S-LAB-NET-20
 2023 09 25 0 8 0
 OPERATOR(S) S.Z.

Net 20 µm

- Decknet
- Deployed at sea

SAMPLE SPLITTING	# of cod-ends [] 1 [] 2	Total volume [] 1600 mL	Aliquots vol. [] 200 mL			
PROTOCOLS	Barcode	Fraction of total volume	Aliquot Volume (mL)	Barcode	Fraction of total volume	Aliquot Volume (mL)
<input checked="" type="radio"/> S20 Cryo-5mL LN2 #1		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL 80 mL <input checked="" type="checkbox"/> 15 mn		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL 90 mL <input checked="" type="checkbox"/> 15 mn
FCAM20 Bottle-250mL LIVE		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL			
E20 Falcon-15mL + 15mL ETOH FRZ -20°C		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL 80 mL			
S20-L Falcon-5mL FRZ -20°C + 5 mL NucleoProtect		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL <input checked="" type="checkbox"/> 15 mn 90 mL			
MB20 Vial-4mL FRZ -20°C		<input type="checkbox"/> 1/8	<input checked="" type="checkbox"/> 200 mL			
FM20 Falcon-50mL FRG +4°C Prealiquoted PFA+GLUTA store at -20°C		<input checked="" type="checkbox"/> 45 mL			<input checked="" type="checkbox"/> 45 mL	



	COMMENTS	COMMENTS
SAMPLE SPLITTING		
PROTOCOLS		
S20 Cryo-5mL LN2 #1		
FCAM20 Bottle-250mL LIVE		
E20 Falcon-15mL FRZ -20°C		
S20-L Falcon-5mL FRZ -20°C		
MB20 Vial-4mL FRZ -20°C		
FM20 Falcon-50mL FRG +4°C		



LOG-SAMPLES_ YYYY MM DD # # # _STATION- # # # _S-LAB-NET-200

OPERATOR(S) S.Z.

**Horizontal
WPII-200**

SAMPLE SPLITTING	COD-END #1					
PROTOCOLS	Barcode	Fraction of total volume	Aliquot Volume (mL)			
F200 Bottle-250mL + borax/formol RT >10°C		<input type="checkbox"/> 1 (100%)	<input type="checkbox"/> 250 mL			
SAMPLE SPLITTING	COD-END #2	Total volume	Aliquots vol.			
PROTOCOLS	Barcode	Fraction of total volume	Aliquot Volume (mL)	Barcode	Fraction of total volume	Aliquot Volume (mL)
S200 Cryo-5mL LN2 #1		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL 90 <input checked="" type="checkbox"/> 15 mn		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL 130ml <input checked="" type="checkbox"/> 15 mn
S200-L Falcon-5mL FRZ -20°C + 5mL Nucleo		<input type="checkbox"/> 1/8	<input type="checkbox"/> 200 mL 140 <input checked="" type="checkbox"/> 15 mn			



	COMMENTS	COMMENTS
SAMPLE SPLITTING		
PROTOCOLS		
F200 Bottle-250mL RT >10°C		
SAMPLE SPLITTING		
PROTOCOLS		
S200 Cryo-5mL LN2 #1		
S200-L Falcon-5mL FRZ -20°C		



YYYY MM DD # # #
 LOG-SAMPLES_ **2023** **09** **25** _STATION- **0** **8** **0** _S-LAB-NET-680
 OPERATOR(S) **S.Z.**

Régent 680

SAMPLE SPLITTING	NET TOW #1			NET TOW #2		
	Total volume [] 1600 mL			Total volume [] 1600 mL		
PROTOCOLS	Barcode	Fraction of total volume	Bottle volume (mL)	Barcode	Fraction of total volume	Aliquot Volume (mL)
F680 Bottle-250mL RT >10°C + Borax/Formol	 112557233	<input type="checkbox"/> 50 % <input type="checkbox"/> 100 %	<input type="checkbox"/> 250 mL 160 mL			
F2000 Bottle-250mL RT >10°C + borax/formol	### EPI F2000	hand-picked #ind=	<input type="checkbox"/> 250 mL			
S680-L Falcon-5mL FRZ -20°C + 5mL Nucleoprotect				 112557234	<input type="checkbox"/> 50 % <input type="checkbox"/> 100 %	15 160 mL 15 15 mn



COMMENTS

COMMENTS

SAMPLE
SPLITTING

PROTOCOLS

F680
Bottle-250mL
RT >10°C

F2000
Bottle-250mL
RT >10°C

S680-L
Falcon-5mL
FRZ -20°C